



iDentalglove™



E-Brochure



Talk to Us

 FLEXIBLE

 SAFE

 COMFORTABLE



@smartgloveintl



info@smartgloveintl.com

Smart Glove International Pte. Ltd



Available in a variety of scents to choose from
Grape, Mint, Bubble Gum, Green Tea, Vanilla and Orange

*All pictures shown are for illustration purpose only. Actual product may vary due to product enhancement.

POLYCHLOROPRENE EXAMINATION GLOVES C-Series

Ultra-Soft Polychloroprene Medical Examination Glove

Material	Polychloroprene
Length	9-inch
Finishing	Powder Free
Design & Features	Finger Textured, Ambidextrous, Full Textured, Beaded Cuff
Colours	Green / Blue / Black / Pink / Brown
Quality Standards	ASTM D6977, EN455 1-4, PPE Cat III EN ISO 374-1 & 5, EN420
Weight (M Size)	3.5g - 6.0g (3mil to 6mil)
Highlights	Provides a very comfortable and soft feel to the wearers with close similarity to natural rubber Low modulus and high elongation help prevent hand fatigue Suitable for dental, medical personnel, veterinary operations, and pharmaceutical

NITRILE EXAMINATION GLOVES S-Series

Super-Soft Nitrile Medical Examination Glove

Material	S1 Nitrile
Length	9-inch / 12-inch
Finishing	Powder Free
Design & Features	Finger Textured, Ambidextrous, Beaded Cuff
Colours	Blue / White / Black / Pink / Orange
Quality Standards	ASTM D6319, EN455 1-4, PPE Cat III EN ISO 374-1 & 5, EN420
Weight (M Size)	3.5g - 8.0g (3mil to 6mil)
Highlights	S-Series offers superior comfort and easy donning to the wearer Offers Accelerator, Sulphur, and Zinc Free option for advanced allergy protection Suitable for medical personnel, veterinary operations, pharmaceutical, and general hygiene use

iDental *glove*™

2018 Asia-Pacific Global Medical Glove Frost & Sullivan's Technology Innovation Award

“ Our award-winning exam gloves are equipped with the latest improved technology, designed to give superior comfort and easy donning. We offer a wide range of exam gloves designed to give superior comfort and maximum barrier protection by exceeding international standards.

With our new, improved technology, we offer ultra-soft polychloroprene copolymer gloves tailored to provide a comfort level similar to natural rubber. Besides providing cross-contamination protection during medical procedures, we cater to personnel who handle chemotherapy drugs as well. We offer a series of scented and unscented gloves crafted with advance formula. Our S-Series and C-series technology produce lighter, softer, and durable gloves. The low modulus and high elongation quality helps prevent hand fatigue for extended wear.



”

FROST & SULLIVAN BEST PRACTICES AWARD



Multiple Use Cases

Suitable for medical personnel, veterinary operations, pharmaceutical & general hygiene use

Allergy-free

Our gloves are made with nitrile - deeming them suitable for users with latex allergies



Comfortable

Providing dexterity and strength, our industry-leading gloves are easy to don and provide superb comfort

Uncompromised quality

Our gloves comply strictly with international quality standards, bringing you assurance and peace of mind

Standards & Compliances

ASTM D6319

Standard Specification for Nitrile Examination Gloves for Medical Application

ASTM D6977

Standard Specification for Polychloroprene Examination Gloves for Medical Application

EN 455

Requirements and Test for Single use medical gloves

Regulation (EU) 2016/425

Personal Protective Equipment - Category III Product

EN 420:2003 + A1

Protective Gloves: General Requirements and Test Methods

ISO 374-1, EN 374-5

Protective Gloves against Dangerous Chemicals and Micro-organisms -Type B

*Other specific local regulatory requirements can be fulfilled upon request.

why iDentalglove™



Superior comfort and easy donning to the wearer



Enhanced performance



Excellent sensitivity and grip



Advanced allergy protection

Applications



Dental



Medical



Pharmaceutical



Veterinary

Our Brands:

GEN-X®

iDentalglove™

SmartMedic™

MetalMan™

iChefglove™

SmartClean™

iSafe™