



SmartMedic™



Our surgical gloves are designed and custom-made to meet the high expectations of surgeons. We offer various types of materials in order to meet the different requirement for different types of surgeries.

➤➤➤ Provides superior softness to aid easy donning and to reduce hand fatigue

➤➤➤ Sterile Exam gloves are available in singles or pairs per pouch packaging

➤➤➤ Tested safe for use with chemotherapy drugs per global Chemotherapy Permeation Test standard



SmartMedic™
Sterile Surgical Glove



Material	Nitrile/ Polychloroprene/ Polyisoprene/ Natural Rubber
Length	12-inch
Finishing	Powdered/ Powder Free
Design & Features	Hand Specific, Micro Textured, Beaded Cuff
Colours	Blue/ White/ Green
Quality Standards	ASTM D3577, ASTM D6978, EN455 1-4, EN420
Weight (M Size)	8.6g - 12.5g (6mil underglove - 8mil)
Highlights	The hand specific gloves are designed considering the anatomy of human hand in order to provide more comfort and less fatigue under prolonged time of wearing Damp Hand Don feature available Among the various surgical gloves, the surgical glove made out of Polyisoprene is the most comfortable to the surgeon

SmartMedic™
Sterile Examination Glove



Material	Nitrile
Length	9-inch, 12-inch
Finishing	Powder Free
Design & Features	Ambidextrous/ Finger Textured or Full Textured / Beaded Cuff
Colours	Blue/ White
Quality Standards	ASTM D6319, ASTM D6978, EN455 1-4, EN420
Weight (M Size)	4.0g - 6.6g (3.5mil - 6mil)



Standards & Compliances

ASTM D3577

Standard Specification for Rubber Surgical Gloves

ASTM D6978-05 (R2019)

Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs

EN 455

Requirements and Test for Single use medical gloves

Regulation (EU) 2016/425

Personal Protective Equipment

EN ISO 21420

Protective Gloves: General Requirements and Test Methods

*Other specific local regulatory requirements can be fulfilled upon request.

Why SmartMedic™



Excellent sensitivity and grip



Superior comfort and easy donning to the wearer



Advanced allergy protection



Protection against contamination

Applications



Surgical Procedure



Cleanroom environment



Controlled environment

Our Brands:

GEN-X®

iDentalglove™

SmartMedic™

MetalMan™

iChefglove™

SmartClean™

iSafe™