



 **iDental** *glove*™



**POLYCHLOROPRENE EXAMINATION GLOVES C-Series** Ultra-Soft Polychloroprene Medical Examination Glove

<b>Material</b>	Polychloroprene
<b>Length</b>	9-inch
<b>Finishing</b>	Powder Free
<b>Design &amp; Features</b>	Finger Textured, Ambidextrous, Full Textured, Beaded Cuff
<b>Colours</b>	Green / Blue / Black / Pink / Brown
<b>Quality Standards</b>	ASTM D6977, EN455 1-4, PPE Cat III EN ISO 374-1 & 5, EN420
<b>Weight (M Size)</b>	3.5g - 6.0g (3mil to 6mil)
<b>Highlights</b>	Provides a very comfortable and soft feel to the wearers with close similarity to natural rubber Low modulus and high elongation help prevent hand fatigue Suitable for dental, medical personnel, veterinary operations, and pharmaceutical

**NITRILE EXAMINATION GLOVES S-Series** Super-Soft Nitrile Medical Examination Glove

<b>Material</b>	S1 Nitrile
<b>Length</b>	9-inch / 12-inch
<b>Finishing</b>	Powder Free
<b>Design &amp; Features</b>	Finger Textured, Ambidextrous, Beaded Cuff
<b>Colours</b>	Blue / White / Black / Pink / Orange
<b>Quality Standards</b>	ASTM D6319, EN455 1-4, PPE Cat III EN ISO 374-1 & 5, EN420, FDA 21 CFR 177.2600
<b>Weight (M Size)</b>	3.5g - 8.0g (3mil to 6mil)
<b>Highlights</b>	S-Series offers superior comfort and easy donning to the wearer Offers Accelerator, Sulphur, and Zinc Free option for advanced allergy protection Suitable for medical personnel, veterinary operations, pharmaceutical, and general hygiene use

# iDentalglove™

## 2018 Asia-Pacific Global Medical Glove Frost & Sullivan's Technology Innovation Award

“ Our award-winning exam gloves are equipped with the latest improved technology, designed to give superior comfort and easy donning. We offer a wide range of exam gloves designed to give superior comfort and maximum barrier protection by exceeding international standards.

With our new, improved technology, we offer ultra-soft polychloroprene copolymer gloves tailored to provide a comfort level similar to natural rubber. Besides providing cross-contamination protection during medical procedures, we cater to personnel who handle chemotherapy drugs as well. We offer a series of scented and unscented gloves crafted with advance formula. Our S-Series and C-series technology produce lighter, softer, and durable gloves. The low modulus and high elongation quality helps prevent hand fatigue for extended wear.



”

### FROST & SULLIVAN BEST PRACTICES AWARD



2018  
GLOBAL MEDICAL GLOVE  
TECHNOLOGY INNOVATION  
AWARD



2019  
ASIA-PACIFIC  
ULTRA-LOW WEIGHT NITRILE GLOVE  
PRODUCT LEADERSHIP AWARD



2020  
ASIA-PACIFIC  
METAL DETECTABLE NITRILE GLOVES  
FOR THE FOOD INDUSTRY  
NEW PRODUCT INNOVATION AWARD



2020  
GLOBAL MEDICAL GLOVES  
TECHNOLOGY INNOVATION  
LEADERSHIP AWARD

#### Multiple Use Cases

Suitable for medical personnel, veterinary operations, pharmaceutical & general hygiene use

#### Allergy-free

Our gloves are made with nitrile - deeming them suitable for users with latex allergies



#### Comfortable

Providing dexterity and strength, our industry-leading gloves are easy to don and provide superb comfort

#### Uncompromised quality

Our gloves comply strictly with international quality standards, bringing you assurance and peace of mind

**Standards & Compliances**

<p><b>ASTM D6319</b></p>	<p><b>ASTM D6977</b></p>	<p><b>EN 455</b></p>	<p><b>Regulation (EU) 2016/425</b></p>
<p>Standard Specification for Nitrile Examination Gloves for Medical Application</p>	<p>Standard Specification for Polychloroprene Examination Gloves for Medical Application</p>	<p>Requirements and Test for Single use medical gloves</p>	<p>Personal Protective Equipment - Category III Product</p>
<p><b>EN 420:2003 + A1</b></p>	<p><b>ISO 374-1, EN 374-5</b></p>	<p><b>Regulation (EU) No 10/2011</b></p>	<p><b>FDA 21 CFR Part 177.2600</b></p>
<p>Protective Gloves: General Requirements and Test Methods</p>	<p>Protective Gloves against Dangerous Chemicals and Micro-organisms -Type B</p>	<p>Plastic Materials and Articles intended to come into contact with food</p>	<p>Substances for Use as Basic Components of Single and Repeated Use Food Contact Surfaces</p>

\*Other specific local regulatory requirements can be fulfilled upon request.

**why iDentalglove™**

			
<p>Superior comfort and easy donning to the wearer</p>	<p>Enhanced performance</p>	<p>Excellent sensitivity and grip</p>	<p>Advanced allergy protection</p>

**Applications**

			
<p>Dental</p>	<p>Medical</p>	<p>Pharmaceutical</p>	<p>Veterinary</p>

**Our Brands:**










Contact us for more information.  
**Smart Glove International Pte. Ltd**

(W) [www.smartgloveintl.com](http://www.smartgloveintl.com)  
✉ [info@smartgloveintl.com](mailto:info@smartgloveintl.com)

(T) +603 5614 7611  
f @smartgloveintl

(F) +603 8601 9488  
t @smartgloveintl