

# Breathfidence™ Face Mask

## ASTM Level-2 [3-ply] Face Mask

### Material & Product Description

Non-Woven PP & Meltblown Filter,  
Layer 1/2/3: 28/22/23 g/m2 each,  
PP Covered Iron Wire as Nose Piece

### Color

Blue and Black

### Product Code

M-000001244-00051 - BL  
F-000002518-00051-BK

### Packing Style

50 pcs per dispenser  
10 dispenser per carton

### Loading Capacity

1,300 Cartons per 1 x 20'GP  
2,600 Cartons per 1 X 40'GP  
3,200 Cartons per 1 X 40'HQ

### Width

86 - 95mm

### Length of earloop

145-185mm

### Product Performance Factors (Nelson Lab Results)

Face Mask Material Requirement by Performance Levels

ASTM Test Standard	Level 1	Level 2	Level 3
BFE (Bacteria Filtration Efficiency) ASTM F2101-07	≥ 95%	≥ 98%	≥ 98%
PFE (Particulate Filtration Efficiency) ASTM F2299	≥ 95%	≥ 98%	≥ 98%
Delta P (Differential Pressure)	< 4.0	< 5.0	< 5.0
Fluid Resistance ASTM F11862 (mm Hg)	80	120	160
Flame Spread - 16 CFR part 1610 (a)	Class 1	Class 1	Class 1



### Dispenser Dimension

188mm x 100mm x 100mm



### Carton Dimension

510mm x 200mm x 210mm

